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Referências

- Adcock DM, Strandberg K, Shima M, Marlar RA. Advantages, disadvantages and optimization of one-stage and chromogenic factor activity assays in haemophilia A and B. *Int J Lab Hematol* 2018; 40(6): 621-629.
- Aznar JA, Cabrera N, Matysiak M, Zawilska K, Gercheva L, Antonov A, Montañés M, Páez AM, Lissitchkov T. Pharmacokinetic study of a high-purity factor IX concentrate (factor IX Grifols) with a 6-month follow up in previously treated patients with severe haemophilia b. *Haemophilia* 2009; 15(6): 1243-1248.
- Baglin T, Huntington JA, Koch A, Mocanu I, Makhaldiani L. Serpin-PC in persons with severe hemophilia (PwH): Updated results from a multicenter multi-part, first-in-human study. *Blood* 2023; 142(Supplement 1) :2619.
- Baker P, Platton S, Gibson C, Gray E, Jennings I, Murphy P, Laffan M. Guidelines on the laboratory aspects of assays used in haemostasis and thrombosis. *Br J Haematol* 2020; 191(3): 347-362.
- Barrowcliffe TW, Raut S, Sands D, Hubbard AR. Coagulation and chromogenic assays of factor VIII activity: General aspects, standardization, and recommendations. *Semin Thromb Hemost* 2002; 28(3): 247-256.

- Bowyer A, Gray E, Lowe A, Murphy P, Platton S, Riddell A, Chowdary P, Lester W, Jenkins PV. Laboratory coagulation tests and recombinant porcine factor VIII: A United Kingdom Haemophilia Centre Doctors' Organisation guideline. *Haemophilia* 2022; 28(3): 515-519.
- Bowyer A, Kitchen S, Maclean R. Effects of emicizumab on APTT, one-stage and chromogenic assays of factor VIII in artificially spiked plasma and in samples from haemophilia A patients with inhibitors. *Haemophilia* 2020; 26(3): 536-542.
- Bowyer A, Kitchen S, Maclean R. Measurement of antifactor VIII antibody titre in the presence of emicizumab; Use of chromogenic Bethesda assays. *Int J Lab Hematol* 2021; 43(4): O204-O206.
- Bowyer AE, Duncan EM, Antovic JP. Role of chromogenic assays in haemophilia A and B diagnosis. *Haemophilia* 2018; 24(4): 578-583.
- Bowyer AE, Ezban M, Kitchen S. Measuring the FVIII mimetic activity of the new bispecific antibody, Mim8, in severe haemophilia A plasma using APTT and one-stage FVIII assays. *Res Pract Thromb Haemost* 2021; 5: Abstract PB0680.
- Bowyer AE, Goodeve AC, Liesner R, Mumford AD, Kitchen S, Makris M. p.Tyr365Cys change in factor VIII: haemophilia A, but not as we know it. *Br J Haem* 2011; 154(5): 618-625.
- Bowyer AE, Gosselin RC. Factor VIII and factor IX activity measurements for hemophilia diagnosis and related treatments. *Semin Thromb Hemost* 2022; 49(06): 609-620.
- Bowyer AE, Hillarp A, Ezban M, Persson P, Kitchen S. Measuring factor IX activity of nonacog beta pegol with commercially available one-stage clotting and chromogenic assay kits: A two-center study. *J Thromb Haemost* 2016; 14(7): 1428-1435.
- Bowyer AE, Kitchen S, Ezban M. The effect of a next generation factor VIII mimetic bispecific antibody (Mim8) on assays of factor VIII activity and thrombin generation. *J Thromb Haemost* 2023; 21(3): 480-487.
- Bowyer AE, Kitchen S, Maclean RM. Effects of emicizumab on APTT, one-stage and chromogenic assays of factor VIII in artificially spiked plasma and in samples from haemophilia A patients with inhibitors. *Haemophilia* 2020; 26(3): 536-542.
- Chhabra ES, Sadeghi-Khomami A, Liu M, Young G, Pipe SW, Ozelo MC, Le-Camus C, Toh M, Lima-Montalvo SA, Demissie M. Global comparative antithrombin (AT) field study: Impact of laboratory assay variability on the assessment of AT activity measurement (abstract PP-072). *Haemophilia* 2024; 30: 3-223.
- Clinical and Laboratory Standards Institute. Determination of coagulation factor activities using the one-stage clotting assay, 2nd edition. CLSI standard H48. 2016. <https://clsi.org/standards/products/hematology/documents/h48/>.
- De Wolf D, Singh K, Chuah MK, VandenDriessche T. Hemophilia gene therapy: The end of the beginning? *Hum Gene Ther* 2023; 34(17-18): 782-792.
- Devreese K, Wijns W, Combes I, Van kerckhoven S, Hoylaerts MF. Thrombin generation in plasma of healthy adults and children: chromogenic versus fluorogenic thrombogram analysis. *Thromb Haemost* 2007; 98(3): 600-613.
- Druzgal CH, Kizilacak H, Brown J, Sennett M, Young G. Neutralizing antidrug antibody to emicizumab in a patient with severe hemophilia A with inhibitors: New case with detailed laboratory evaluation. *J Thromb Haemost* 2020; 18(9): 2205-2208.
- European Medicines Agency (EMA). Hemlibra (emicizumab) Summary of Product Characteristics. 2018. https://www.ema.europa.eu/en/documents/product-information/hemlibra-epar-product-information_en.pdf.
- European Medicines Agency (EMA). NovoEight (turoctocog alfa) Summary of Product Characteristics. 2021. https://www.ema.europa.eu/en/documents/product-information/novoeight-epar-product-information_en.pdf.
- European Medicines Agency (EMA). Hemgenix (etranacogene dezaparvovec) Summary of Product Characteristics. 2024. https://www.ema.europa.eu/en/documents/product-information/hemgenix-epar-product-information_en.pdf.
- European Medicines Agency (EMA). Nuwiq (simoctocog alfa) Summary of Product Characteristics. 2022. https://www.ema.europa.eu/en/documents/product-information/nuwiq-epar-product-information_en.pdf.
- Foley JH, Shehu E, Riddell A, Gray E, Goodale A, Yu IM et al. Differences in wild-type- and R338L-tenase complex formation are at the root of R338L-factor IX assay discrepancies. *Blood Adv* 2023; 7(3): 458-467.
- Food and Drug Administration (FDA). Xyntha (antihemophilic factor [recombinant]) Prescribing Information. 2020. <https://www.fda.gov/media/70399/download>.

Gray E, Kitchen S, Bowyer AE, Chowdary P, Jenkins PV, Murphy P et al. Laboratory measurement of factor replacement therapies in the treatment of congenital haemophilia: A United Kingdom Haemophilia Centre Doctors' Organisation guideline. *Haemophilia* 2020; 26(1): 6-16.

Gritsch H, Romeda-Finger S, Scheiflinger F, Turecek PL. Potency assignment and measurement of recombinant FIX activity in human plasma – impact of aPTT reagents on the 1-stage clotting assay. *Haemophilia* 2014; 20(s3): 37 (abstract).

Ingerslev J, Jankowski MA, Weston SB, Charles LA. Collaborative field study on the utility of a BDD factor VIII concentrate standard in the estimation of BDDr factor VIII:C activity in hemophilic plasma using one-stage clotting assays. *J Thromb Haemost* 2004; 2(4): 623-628.

Jenkins PV, Bowyer AE, Burgess C, Gray E, Kitchen S, Murphy P et al. Laboratory coagulation tests and emicizumab treatment A United Kingdom Haemophilia Centre Doctors' Organisation guideline. *Haemophilia* 2020; 26: 151-155.

Jennings I, Kitchen DP, Woods TA, Kitchen S, Walker ID. Emerging technologies and quality assurance: the United Kingdom National External Quality Assessment Scheme perspective. *Semin Thromb Hemost* 2007; 33(3): 243-249.

Jiang M, Yang F, Jiang Y, Cheng L, Han J, Yi J et al. Safety and efficacy of an anti-human APC antibody for prophylaxis of congenital factor deficiencies in preclinical models. *Blood* 2023; 142(12): 1071-1081.

Kershaw GW, Dissanayake K, Chen VM, Khoo T. Evaluation of chromogenic FIX assays by automated protocols. *Haemophilia* 2018; 24(3): 492-501.

Kitchen S, Beckmann H, Katterle Y, Bruns S, Tseneklidou-Stoeter D, Maas Enriquez M. BAY 81-8973, a full-length recombinant factor VIII: results from an International comparative laboratory field study. *Haemophilia* 2016; 22(3): e192-e199.

Kitchen S, Jennings I, Makris M, Kitchen DP, Woods TAL, Walker ID. Clotting and chromogenic factor VIII assay variability in post-infusion and spiked samples containing full-length recombinant FVIII or recombinant factor VIII Fc fusion protein (rFVIII-Fc). *Int J Lab Hematol* 2019; 41(2): 176-183.

Kitchen S, Jennings I, Makris M, Kitchen DP, Woods TAL, Walker ID. Factor VIII assay variability in postinfusion samples containing full length and B-domain deleted FVIII. *Haemophilia* 2016; 22(5): 806-812.

Kitchen S, Katterle Y, Beckmann H, Maas Enriquez M. Chromogenic assay for BAY 81-8973 potency assignment has no impact on clinical outcome or monitoring in patient samples. *J Thromb Haemost* 2016; 14(6): 1192-1199.

Kitchen S, Kershaw GW, Tiefenbacher S. Recombinant to modified factor VIII and factor IX – chromogenic and one-stage assay issues. *Haemophilia* 2016; 22: 72-77.

Klukowska A, Szczepański T, Vdovin V, Knaub S, Jansen M, Liesner R. Novel, human cell line-derived recombinant factor VIII (Human-cl rhFVIII, Nuwiq®) in children with severe haemophilia A: efficacy, safety and pharmacokinetics. *Haemophilia* 2016; 22(2): 232-239.

Lenting PJ, Denis CV, Christophe OD. Emicizumab, a bispecific antibody recognising coagulation factors IX and X: how does it actually compare to factor VIII? *Blood* 2017; 130: 2463-2468.

Lissitchkov T, Hampton K, von Depka M, Hay C, Rangarajan S, Tuddenham E et al. Novel, human cell line-derived recombinant factor VIII (human-cl rhFVIII; Nuwiq®) in adults with severe haemophilia A: efficacy and safety. *Haemophilia* 2016; 22(2): 225-231.

Mahlangu J, Kaczmarek R, von Drygalski A, Shapiro S, Chou S, Ozelo MC et al. Two-year outcomes of valoctocogene roxaparvovec therapy for hemophilia A. *New Eng J Med* 2023; 388(8): 694-705.

Mahlangu JN, Ahuja SP, Windyga J, Church N, Shah A, Schwartz L. BAY 81-8973, a full-length recombinant factor VIII for the treatment of hemophilia A: product review. *Ther Adv Hematol* 2018; 9(7): 191-205.

Mast AE, Ruf W. Regulation of coagulation by tissue factor pathway inhibitor: Implications for hemophilia therapy. *J Thromb Haemost* 2022; 20(6): 1290-1300.

Matino D, Acharya S, Palladino A, Hwang E, McDonald R, Taylor CT et al. Efficacy and safety of the anti-tissue factor pathway inhibitor marstacimab in participants with severe hemophilia without inhibitors: Results from the phase 3 Basis trial. *Blood* 2023; 142(Supplement 1): 285.

- Miesbach W, Meijer K, Coppens M, Kampmann P, Klamroth R, Schutgens R et al. Gene therapy with adeno-associated virus vector 5-human factor IX in adults with hemophilia B. *Blood* 2018;131(9):1022-31.
- Miller CH, Boylan B, Payne AB, Driggers J, Bean CJ. Validation of the chromogenic Bethesda assay for factor VIII inhibitors in hemophilia A patients receiving Emicizumab. *Int J Lab Hematol* 2021; 43(2): e84-e86.
- Miller CH, Rice AS, Boylan B, Shapiro AD, Lentz SR, Wicklund BM et al. Comparison of clot-based, chromogenic, and fluorescence assays for measurement of factor VIII inhibitors in the U.S. Hemophilia Inhibitor Research Study. *J Thromb Haemost* 2013; 11(7): 1300-1309.
- Morfini M, Cinotti S, Bellatreccia A, Paladino E, Gringeri A, Mannucci PM. A multicenter pharmacokinetic study of the B-domain deleted recombinant factor VIII concentrate using different assays and standards. *J Thromb Haemost* 2003; 1(11): 2283-2289.
- Moser KA, Smock KJ. Direct oral anticoagulant (DOAC) interference in hemostasis assays. *Hematology Am Soc Hematol Educ Program* 2021; 2021(1): 129-133.
- Nederlof A, Kitchen S, Meijer P, Cnossen MH, Ali Pour N, Kershaw GW et al. Performance of factor IX extended half-life product measurements in external quality control assessment programs. *J Thromb Haemost* 2020; 18: 1874-1883.
- Ninivaggi M, de Laat-Kremers R, Tripodi A, Wahl D, Zuily S, Dargaud Y et al. Recommendations for the measurement of thrombin generation: Communication from the ISTH SSC Subcommittee on Lupus Anticoagulant/Antiphospholipid Antibodies. *J Thromb Haemost* 2021; 19(5): 1372-1378.
- Nogami K, Shima M. Current and future therapies for haemophilia—Beyond factor replacement therapies. *Br J Haematol* 2023; 200(1): 23-34.
- Pasi KJ, Lissitchkov T, Mamonov V, Mant T, Timofeeva M, Bagot C et al. Targeting of antithrombin in hemophilia A or B with investigational siRNA therapeutic fitusiran-Results of the phase 1 inhibitor cohort. *J Thromb Haemost* 2021; 19(6): 1436-1446.
- Pavlova A, Delev D, Pezeshkpoor B, Muller J, Oldenburg J. Haemophilia A mutations in patients with non-severe phenotype associated with a discrepancy between one-stage and chromogenic factor VIII activity assays. *Thromb Haemost* 2014; 111(5): 851-861.
- Peyvandi F, Oldenburg J, Friedman KD. A critical appraisal of one-stage and chromogenic assays of factor VIII activity. *J Thromb Haemost* 2016; 14(2): 248-261.
- Pipe SW, Leebeek FW, Recht M, Key NS, Castaman G, Miesbach W et al. Gene therapy with etranacogene dezaparvovec for hemophilia B. *New Engl J Med* 2023; 388(8): 706-718.
- Pipe SW, Sadeghi-Khomami A, Konkle BA, Kitchen S, Negrier C, Liu M, et al. A global comparative field study to evaluate the factor VIII activity of efanesoctocog alfa by one-stage clotting and chromogenic substrate assays at clinical haemostasis laboratories. *Haemophilia*. 2024;30(1):214-23.
- Pittman DD, Carrier C, Soares H, McKay J, Tan CY, Liang JZ et al. Field study and correlative studies of factor IX variant FIX-R338L in participants treated with fidanacogene elaparvovec. *Thromb Haemost* 2024; 124(10): 912-921.
- Platton S, Raheja P, Dale C, Guy S, Yartey N, Bowyer A. Evaluation of one-stage and chromogenic assays for the laboratory measurement of factor VIII activity following valoctocogene roxaparvovec infusion. *Haemophilia* 2024; 30(5): 1221-1224.
- Pouplard C, Trossaert M, A LEQ, Delahousse B, Giraudeau B, Gruel Y. Influence of source of phospholipids for APTT-based factor IX assays and potential consequences for the diagnosis of mild haemophilia B. *Haemophilia* 2009; 15(1): 365-368.
- Prince RE, Schaeper U, Dames S, Calzavarini S, Quarroz C, Reina Caro MD et al. Targeting protein S using small interfering RNA is well tolerated and protects mice with hemophilia A from acute hemarthrosis. *Blood* 2020; 136(Supplement 1): 20-21.
- Rangarajan S, Walsh L, Lester W, Perry DJ, Madan B, Laffan M. AAV5-factor VIII gene transfer in severe hemophilia A. *N Engl J Med* 2017; 377(26): 2519-2530.
- Rezaie AR, Giri H. Anticoagulant and signaling functions of antithrombin. *J Thromb Haemost* 2020; 18(12): 3142-3153.

Robinson M, George L, Carr ME, Samuelson-Jones BJ, Arruda VR, Murphy JE et al. Factor IX assay discrepancies in the setting of liver gene therapy using a hyperfunctional variant factor IX-Padua. *J Thromb Haemost* 2021; 19(5):1212-1218.

Rosen S, Tiefenbacher S, Robinson M, Huang M, Srimani J, Mackenzie D et al. Activity of transgene-produced B-domain-deleted factor VIII in human plasma following AAV5 gene therapy. *Blood* 2020; 136(22): 2524-2534.

Rosen S. Assay of factor VIII:C with a chromogenic substrate. *Scand J Haematol* 1984; 33(suppl 40): 139-145.

Ruinemans-Koerts J, Peterse-Stienissen I, Verbruggen B. Non-parallelism in the one-stage coagulation factor assay is a phenomenon of lupus anticoagulants and not of individual factor inhibitors. *Thromb Haemost* 2010; 104(5): 1080-1082.

Simioni P, Tormene D, Tognin G, Gavasso S, Bulato C, Iacobelli NP et al. X-linked thrombophilia with a mutant factor IX (factor IX Padua). *N Engl J Med* 2009; 361: 1671-1675.

Sommer JM, Buyue Y, Bardan S, Peters RT, Jiang H, Kamphaus GD et al. Comparative field study: impact of laboratory assay variability on the assessment of recombinant factor IX Fc fusion protein (rFIXFc) activity. *Thromb Haemost* 2014; 112(5): 932-940.

Srivastava A, Santagostino E, Dougall A, Kitchen S, Sutherland M, Pipe SW et al. WFH Guidelines for the Management of Hemophilia, 3rd edition. *Haemophilia* 2020; 26(suppl 6): 1-158.

Tiefenbacher S, Clausen WHO, Hansen M, Luthoft R, Ezban M. A field study evaluating the activity of N8-GP in spiked plasma samples at clinical haemostasis laboratories. *Haemophilia* 2019; 25: 893-901.

Truedsson Å, Schmidt DE, Strålfors A, Soutari N, Norberg E, Letelier A et al. One-stage versus chromogenic factor IX activity in haemophilia B [abstract]. *Res Pract Thromb Haemost* 2020;4.

Valsecchi C, Gobbi M, Beeg M, Adams T, Castaman G, Schiavone L et al. Characterization of the neutralizing anti-emizumab antibody in a patient with hemophilia A and inhibitor. *J Thromb Haemost* 2021; 19(3): 711-718.

Van Cott EM, Khor B, Zehnder JL. Factor V Leiden. *Am J Hematol* 2016; 91(1): 46-49.

Verbruggen B, Giles AR, Samis J, Verbeek K, Mensink E, Novakova I. The type of factor VIII deficient plasma used influences the performance of the Nijmegen modification of the Bethesda assay for factor VIII inhibitors. *Thromb Haemost* 2001; 86: 1435-1439.

Verbruggen B, Novakova I, Wessels H, Boezeman J, van den Berg M, Mauser-Bunschoten E. The Nijmegen modification of the Bethesda assay for factor VIII:C inhibitors: improved specificity and reliability. *Thromb Haemost* 1995; 72(2): 247-251.

Viuff D, Barrowcliffe TW, Saugstrup T, Ezban M, Lillicrap D. Viuff D, Barrowcliffe T, Saugstrup T, Ezban M, Lillicrap D. International comparative field study of N8 evaluating factor VIII assay performance. *Haemophilia* 2011; 17(4): 695-702.

Wilmot HV, Hogwood J, Gray E. Recombinant factor IX: discrepancies between one-stage clotting and chromogenic assays. *Haemophilia* 2014; 20(6): 981-987.

Wilson HP, Pierre A, Paysse AL, Kumar N, Cooley BC, Rudra P et al. Protein S antibody as an adjunct therapy for hemophilia B. *Blood Advances* 2024; 8(2): 441-452.

Woodhams B, Girardot O, Blanco MJ, Colesse G, Gourmelin Y. Stability of coagulation proteins in frozen plasma. *Blood Coagul Fibrinolysis* 2001; 12(4): 229-236.

Young G, Lenting PJ, Croteau SE, Nolan B, Srivastava A. Antithrombin lowering in hemophilia: a closer look at fitusiran. *Res Pract Thromb Haemost* 2023; 7(4): 100179.

Zhao Y, Feng G, Feng L. Effects of pre-analytical storage time, temperature, and freeze-thaw times on coagulation factors activities in citrate-anticoagulated plasma. *Ann Transl Med* 2018; 6(23): 456.